

Sizing Guide

Mesial-Distal dimensions in millimeters (mm)



CENTRAL	
Size	Maxillary
0	6.12
1	6.48
2	6.87
3	7.22
4	7.58
5	7.96
6	8.35

LATERAL	
Size	Maxillary
0	4.73
1	5.06
2	5.41
3	5.79
4	6.20
5	6.63
6	6.96

UNIVERSAL	
Size	Mandibular
1	3.80
2	4.09
3	4.40
4	4.73

CUSPID		
Size	Maxillary	Mandibular
0	6.20	5.23
1	6.51	5.50
2	6.84	5.78
3	7.18	6.06
4	7.54	6.37
5	7.91	6.69
6	8.31	7.02

1ST MOLAR		
Size	Maxillary	Mandibular
1	6.52	7.26
2	6.85	7.70
3	7.26	8.16
4	7.62	8.65
5	8.00	9.08
6	8.41	9.54
7	8.83	10.01

1ST MOLAR (N)		
Size	Maxillary	Mandibular
1	6.04	6.85
2	6.35	7.20
3	6.73	7.63
4	7.07	8.09
5	7.42	8.49
6	7.79	8.92
7	8.18	9.36

2ND MOLAR		
Size	Maxillary	Mandibular
1	8.49	9.04
2	9.00	9.50
3	9.45	9.97
4	9.92	10.47
5	10.42	10.89
6	10.83	11.33
7	11.27	11.78

2ND MOLAR (N)		
Size	Maxillary	Mandibular
1	7.63	8.17
2	8.09	8.58
3	8.49	9.03
4	8.92	9.47
5	9.37	9.85
6	9.74	10.25
7	10.13	10.66

PERMANENT MOLAR		
Size	Maxillary	Mandibular
2	10.42	10.85
3	10.85	11.30
4	11.23	11.70
5	11.62	12.10
6	12.03	12.53
7	12.45	12.97